

PRIVILEGE LIST - AUDIOLOGY

INSTRUCTIONS

Part I, List of Privileges: The **PROVIDER** enters the appropriate code number in the block marked **REQUESTED** for each privilege. Each block must have a code number. The provider signs and dates the form and sends it to his or her supervisor.

Part II, Supervisor's Recommendation: The supervisor reviews the requested privileges and checks one of the blocks, signs and dates the form and sends it to the Credentials Function Chairperson.

The CREDENTIALS FUNCTION CHAIRPERSON: Enters the appropriate code number in the block marked **APPROVED** for each privilege. Each block must have a code number.

The MEDICAL FACILITY COMMANDER (MFC): Acting on the recommendations of the credentials function, awards, renews, and revises individual provider privileges.

Revisions and Corrections: To change a code number, the originator of the revision forwards a written request to the Credentials Function Chairperson. The function sends its recommendation to the MFC in turn. After approval the Credentials Function Chairperson makes the change in pen, initials it, and notifies the provider in writing of the decision. Any corrections must be initialed in pen. The credentials function may require the provider to complete a new privilege list

- CODES:**
- | | |
|--|---|
| <p>1. Perform unsupervised</p> <p>2. Perform with supervision.</p> | <p>3. Not requested/approved due to lack of facility support</p> <p>4. Not requested and/approved due to lack of expertise.</p> |
|--|---|

NAME OF PRACTITIONER (Last, First, Middle Initial)

NAME OF MEDICAL FACILITY

I. LIST OF PRIVILEGES - AUDIOLOGIST

REQUEST	APPROVED		REQUEST	APPROVED	
		A. BASIC AUDIOMETRY			6. Bone Conduction Implantable hearing aids
		1. Pure Tone and speech audiometry			7. Rehabilitation counseling and training
		2. Acoustic immittance			D. SUPPLEMENTAL PRIVILEGES
		B. ADVANCED AUDIOMETRY			1. Electroneuronography
		1. Audiologic site of lesion battery			2. Cerumen Management
		2. Fitness for duty assessment			3. Otolaryngology resident training
		3. Functional hearing loss assessment			4. Audiometric technician certification course director
		4. Electrocochleography (ECOG)			5. Educational lectures
		5. Audiometric Brainstem Response (ABR)			E. ADDITIONAL PRIVILEGES
		6. Middle Latency Response (MLR)			1. Refer patients to specialists and other providers as warranted by their condition
		7. Long Latency Responses (LLR)			2. Request appropriate laboratory studies (to be interpreted by a physician)
		8. Cognitive Evoked Responses (P-300)			F. OTHER PRIVILEGES (LIST)
		9. Cochlear implant evaluations			
		10. Central auditory assessment			
		11. Electronystagmography			
		12. Posturography			
		13. Rotational chair testing			
		14. Hearing Conservation Program Management			
		15. Otoacoustic emissions testing			
		16. Neonatal screening (Term Nursery & NICU)			
		17. Intraoperative monitoring of cranial nerves			
		C. AMPLIFICATION AND AURAL REHABILITATION			
		1. Hearing aid candidacy assessment, selection, fitting and counseling			
		2. Earmold fabrication and modification			
		3. Electroacoustic measurement of hearing aids			
		4. Real ear measurements			
		5. Programmable hearing aids			

SIGNATURE OF PRACTITIONER

DATE

**Download any U.S. FedForm (free, fillable, savable in Adobe Reader)!
Start with the "Flash Demo" at the top of the following page:
www.usa-federal-forms.com**

**Convert any fillable PDF form to savable (locally, in Adobe Reader):
www.savePDF.com**

**Convert any document (in any format) to PDF fillable and savable:
www.FillinDocs.com**

**All (10's of 1,000's) U.S. Federal Forms already fillable, savable:
www.usa-federal-forms.com**

About the ITAOP/savePDF Method

The traditional Field-by-Field creation process is extremely ineffective and slow.

The only realistic option to create high-quality forms is the Insert-Text-Anywhere-on-Page (ITAOP) method.

The field creation process is about 10,000 times faster than the traditional method; the list of ITAOP features is not even available for the traditional method.

ITAOP savePDF method proved to be very simple and completely reliable for millions of users all over the world (incl. individuals, companies, organizations, government employees).

II. SUPERVISOR'S RECOMMENDATION		
<input type="checkbox"/> RECOMMEND APPROVAL	<input type="checkbox"/> RECOMMEND APPROVAL WITH MODIFICATION <i>(Specify below)</i>	<input type="checkbox"/> RECOMMEND DISAPPROVAL <i>(Specify below)</i>
SIGNATURE OF SUPERVISOR		DATE