

WEEKLY AMBULANCE CHECKLIST

JUMP BAG

<i>IV starting kit</i>	<i>4 x 4's</i>
<i>Pocket mask</i>	<i>ABD pads or field dressings</i>
<i>Kerlex</i>	<i>Cravats</i>
<i>ACE bandages</i>	<i>IV solution (2 x liters NS) (expiration)</i>
<i>Triage cards</i>	<i>Nutrach kit</i>
<i>Tape</i>	<i>Scissors</i>
<i>Stethoscope</i>	<i>Penlight</i>
<i>Oral airways</i>	<i>Latex gloves</i>

O2 BAG

<i>ET tubes</i>	<i>Laryngoscope Handle (function)</i>
<i>Oral airways</i>	<i>Laryngoscope Blades</i>
<i>Nasal airways</i>	<i>ET tube stylettes</i>
<i>O2 masks (adult/peds)</i>	<i>Tape</i>
<i>10 cc syringes</i>	<i>Ambu bag (adult/ped)</i>
<i>Stethoscope</i>	<i>Aviators hose connection</i>
<i>O2 tubing</i>	<i>Portable O2 tank (condition/pressure)</i>

DEFIBRILLATOR

AED

MISCELLANEOUS

<i>Expiration of electrode pads</i>	<i>Spare battery</i>	<i>Gurney function</i>
<i>Paper</i>	<i>Audio cassette</i>	<i>Cervical collars in good repair</i>
<i>Cables (condition)</i>	<i>Memory module</i>	<i>Empty trash bin</i>
<i>Extra electrode pads</i>	<i>Cables (condition)</i>	<i>Backboard set up w/straps & head restraints</i>
<i>Function check</i>	<i>Set of electrodes</i>	
	<i>Function check</i>	

IS THIS AMBULANCE READY FOR ANY CONTINGENCY?

SIGNATURE

DATE

**Download any U.S. FedForm (free, fillable, savable in Adobe Reader)!
Start with the "Flash Demo" at the top of the following page:
www.usa-federal-forms.com**

**Convert any fillable PDF form to savable (locally, in Adobe Reader):
www.savePDF.com**

**Convert any document (in any format) to PDF fillable and savable:
www.FillinDocs.com**

**All (10's of 1,000's) U.S. Federal Forms already fillable, savable:
www.usa-federal-forms.com**

About the ITAOP/savePDF Method

The traditional Field-by-Field creation process is extremely ineffective and slow.

The only realistic option to create high-quality forms is the Insert-Text-Anywhere-on-Page (ITAOP) method.

The field creation process is about 10,000 times faster than the traditional method; the list of ITAOP features is not even available for the traditional method.

ITAOP savePDF method proved to be very simple and completely reliable for millions of users all over the world (incl. individuals, companies, organizations, government employees).