

PATIENT IDENTIFICATION				REQUEST LOCATION	UCA CODE	PROVIDER NUMBER	STAT HEM/UA/COAG				
FMP-SSN				<b>INSTRUCTIONS:</b> 1. NAME FMP AND SSN OF PATIENT MUST BE ON THE LAB SLIP. 2. AT LEAST ONE TEST MUST BE MARKED. 3. ONLY ONE SPECIMEN SOURCE MAY BE MARKED. 4. ONLY ONE COLLECTION MAY BE MARKED. 5. TIME AND DATE MUST BE ENTERED UPON SPECIMEN COLLECTION	PROVIDER STAMP						
NAME											
DOB		SEX	PHONE								
SPONSER		RELATION	SVC					GRADE			
TIME		DATE		LOCATION OF RECORDS							
STAT HEMATOLOGY			STAT COAGULATION		STAT URINALYSIS			COLLECTED BY			
CBC & AUTOMATED DIFFERENTIAL			PT WITH INR		STAT URINALYSIS			ALL AVAILABLE STAT TESTS IN THIS CATEGORY ARE LISTED. ANY EXCEPTIION MUST BE ARRANGED BY CONTACTING THE ON CALL PATHOLOGISTS  CARE PROVIDER COMMENTS			
HCT & HGB ONLY			APTT		STAT URINE PREG TEST						
MANUAL DIFFERENTIAL			FIBRINOGEN		STAT SERUM PREG TEST (Qualitative)						
BODY FLUID EXAMS			THROMBIN TIME		COMMENTS:					1. WARD CLINIC PERSONNEL	
CSF			D-DIMER							2. LAB PERSON IN LAB	
JOINT SYNOVIAL (Include s Crystal Exam)			FDP/FSP					3. WARD ROUNDS			
PERITONEAL			1 : 1 MIX					4. LAB TECH AT WARD/CLINIC			
PLEURAL			TEG					5. OTHER			
PERICARDIAL			OTHER (Specify in Comments)								

WHMC FORM 3265, AUG 98 (IMT-V1) PREVIOUS EDITION WILL BE USED

(1ST COPY - HEME)  
(2ND COPY - COAG)

**Download any U.S. FedForm (free, fillable, savable in Adobe Reader)!  
Start with the "Flash Demo" at the top of the following page:  
[www.usa-federal-forms.com](http://www.usa-federal-forms.com)**

**Convert any fillable PDF form to savable (locally, in Adobe Reader):  
[www.savePDF.com](http://www.savePDF.com)**

**Convert any document (in any format) to PDF fillable and savable:  
[www.FillinDocs.com](http://www.FillinDocs.com)**

**All (10's of 1,000's) U.S. Federal Forms already fillable, savable:  
[www.usa-federal-forms.com](http://www.usa-federal-forms.com)**

---

### **About the ITAOP/savePDF Method**

The traditional Field-by-Field creation process is extremely ineffective and slow.

The only realistic option to create high-quality forms is the Insert-Text-Anywhere-on-Page (ITAOP) method.

The field creation process is about 10,000 times faster than the traditional method; the list of ITAOP features is not even available for the traditional method.

ITAOP savePDF method proved to be very simple and completely reliable for millions of users all over the world (incl. individuals, companies, organizations, government employees).