

(This form is affected by the Privacy Act of 1974- Use blanket PAS DD Form 2005)

PATIENT IDENTIFICATION				LOCATION	UCA CODE	PROVIDER NO	SEROLOGY	
FMP - SSN				INSTRUCTIONS: 1. Name, FMP and SSN of patient must be on the lab slip. 2. Patient ID, Location, UCA code and provider number must be completed. 3. At least one "TEST" must be marked. 4. Only one "SPECIMEN SOURCE" must be marked. 5. Only one "COLLECTED BY" may be marked. 6. Time and date must be entered upon specimen collection.	PROVIDER STAMP			
NAME								
DOB SEX PHONE								
SPONSOR'S RELATION SVC GRADE								
TIME	DATE	LOCATION OF RECORDS						
TEST	RESULTS	TEST	RESULTS	TEST	RESULTS	COLLECTED BY		
Anti-Nuclear Antibody		CMV: <input type="checkbox"/> G <input type="checkbox"/> M		HEP B core IgM AB		1 WARD CLINIC PERSONNEL		
Anti-dsDNA		EBV: <input type="checkbox"/> G <input type="checkbox"/> M		HEP B surface AB		2 LAB PERSON IN LAB		
Anti-ENA PANEL		Infectious MONO		Hepatitis A IgM AB		3 WARD ROUNDS		
Anti-Smith		Rheumatoid factor		Hepatitis C AB		4 LAB PERSON AT		
Anti - RNP		Rubella		ACUTE HEPATITIS PANEL (HBsAG, HBcIgM, HAV, HCV)		OTHER (Specify in comments)		
CYTOPLASTIC AUTO AB PANEL		Rubeola		CHRONIC HEPATITIS PANEL (HBsAG, HBc Total, HCV)		CARE PROVIDER		
Mitochondrial AB		RPR screen		HAPC PANEL				
Smooth muscle AB		Sero TP-PA (MHATP)		H. Pylori AB				
Parietal cell AB		Streptozyme		Other:				
ANCA PANEL		THYROID ANTIBODY PANEL						
ANTI-MPO		Anti - microsomal AB						
ANTI-PR3		Anti - thyroglobulin AB						
ANTICARDIOLIPIN PANEL		Diagnostic Hepatitis						
IgA		HEP B Surface AG						
IgM		HEP B Total core AB						
IgG								

WHMC FORM 3600, 20010201 (IMT-V1) (MTLLI)

**Download any U.S. FedForm (free, fillable, savable in Adobe Reader)!
Start with the "Flash Demo" at the top of the following page:
www.usa-federal-forms.com**

**Convert any fillable PDF form to savable (locally, in Adobe Reader):
www.savePDF.com**

**Convert any document (in any format) to PDF fillable and savable:
www.FillinDocs.com**

**All (10's of 1,000's) U.S. Federal Forms already fillable, savable:
www.usa-federal-forms.com**

About the ITAOP/savePDF Method

The traditional Field-by-Field creation process is extremely ineffective and slow.

The only realistic option to create high-quality forms is the Insert-Text-Anywhere-on-Page (ITAOP) method.

The field creation process is about 10,000 times faster than the traditional method; the list of ITAOP features is not even available for the traditional method.

ITAOP savePDF method proved to be very simple and completely reliable for millions of users all over the world (incl. individuals, companies, organizations, government employees).